

Every year there are about 30,000 cardiac arrests in the UK that happen outside of a hospital environment. For every second that defibrillation is delayed, a casualty's chances of survival falls by about 10%. We've put together a few tips to help you ensure your defibrillator is always ready for use.

4	Set up a regular inspection procedure for your
	defibrillator.

Record your inspections. (We've put together a simple sheet to help you).

Make a note of the expiration date of pads and batteries.

Reorder batteries and pads before they expire.

Follow the guidance in your product manual when cleaning your defibrillator.

Defibrillator ready to use? Y/N	Date	Initial

**Batteries expiration date** 

Pads expiration date

